



Progesterone Rabbit anti Human (AD-PROG3CMO)

Lot: M106240776

Immunogen: 3CMO-BSA

Storage condition: -20°C

Cross reaction:

Liquid Phase: RIA

Working dilution: used 1/75.000

Sensitivity of the standard curve: 4 – 8 pg.

11 α -Hydroxyprogesterone	1.39%
20 β Hydroxyprogesterone	0.4%
17 α -Hydroxyprogesterone	1.12%
Testosterone	0.03%
17 β -oestradiol	less than 0.009%
Estradiol	less than 0.009%
5 β Pregnan 3,20 dione	2.7%
20 Hydroxy delta 4 Pregnen 3,one	4.1%
Pregnelone	3%
Cortisol	less than 0.009%
11-Deoxycorticosterone	2.45%
11-Deoxycortisol	0.12%
Androstenedione	0.13%
DHEA	less than 0.009%
DHEA-SO ₄	less than 0.009%

ALPHADIA S.A/N.V.

Siège d'exploitation : Rue du bosquet, 15A • BE-1435 Mont-Saint-Guibert • BELGIUM

Siège social : Rue René Descartes, 19 • BE-7000 Mons • BELGIUM

RPM Mons • TVA BE 0454.402.636

Tel : +32 (0) 10 68 56 10 • E-mail : contact@alphadia.be • www.alphadia.be